



Broadcast
Bionics

Bionic Talkshow Keyboard Shortcuts

<Duncan Smith>



Bionics Corporation Ltd (trading as Broadcast Bionics),
The Barn, Hurstwood Grange, Hurstwood Lane, Haywards Heath RH17 7QX
T: +44 (0)1444 473 999 F: +44 (0)1444 473 888 E: sales@bionics.co.uk

Bionic Talkshow Keyboard Short cuts

F5 > Puts a call back on hold
F6 > Hangs call up
F7 > Flags a call
F8 > Places the call in the On Air queue
SHIFT + CTRL + F1 > Terminate call on hybrid 1
SHIFT + CTRL + F2 > Terminate call on hybrid 2
SHIFT + CTRL + F3 > Terminate call on hybrid 3
SHIFT + CTRL + F4 > Terminate call on hybrid 4
SHIFT + CTRL + F5 > Terminate call on handset
SHIFT + F1 > Route call from handset to hybrid 1
SHIFT + F2 > Route call from handset to hybrid 2
SHIFT + F3 > Route call from handset to hybrid 3
SHIFT + F4 > Route call from handset to hybrid 4
CTRL F1 > Route from hybrid 1 to handset
CTRL F2 > Route from hybrid 2 to handset
CTRL F3 > Route from hybrid 3 to handset
CTRL F4 > Route from hybrid 4 to handset
CTRL + HOME > Shows
CTRL M to this for create new internal note
Screener > CTRL + 1,2,3,4,5 > Ratings
Screener > CTRL + F, M, O > Gender
Screener > F9 > Dispo page > a-t selection
Screener > F10 > Age page > a-f selection
Screener > F11 > Put call On Air
Screener > F12 > Print function (not just screener?)
Screener > ALT + 0 > 9 tag picker, 0=no tag
All views with Chat > F7 focus in chat window
Screener > F5 > Hang up handset
Screener > F6 > Reject

